

For Research Use Only

# EXD2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85832-3-PBS



## Basic Information

Catalog Number:

85832-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14090

GenBank Accession Number:

BC001962

GeneID (NCBI):

55218

UNIPROT ID:

Q9NVH0

Full Name:

exonuclease 3'-5' domain containing

2

Calculated MW:

70 kD aa, 56 kDa

Purification Method:

Protein A purification

CloneNo.:

250133F10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02142-1, EXD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85832-3-PBS. Detection antibody: 85832-2-PBS. Standard: Ag14090. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02142-2, EXD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85832-3-PBS. Detection antibody: 85832-1-PBS. Standard: Ag14090. Range: 0.781-100 ng/mL.