

For Research Use Only

# Nectin 3 Recombinant antibody

Catalog Number: 85833-4-RR



## Basic Information

Catalog Number:

85833-4-RR

Concentration:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2657

GenBank Accession Number:

NM\_015480.3

GeneID (NCBI):

25945

UNIPROT ID:

Q9NQS3-1

Full Name:

poliovirus receptor-related 3

Calculated MW:

61 kDa

Observed MW:

70-80 kDa

Purification Method:

Protein A purification

CloneNo.:

242804D5

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HCT 116 cells, mouse testis tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85833-4-RR (PVRL3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.