

For Research Use Only

PNO1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85834-2-PBS



Basic Information

Catalog Number:

85834-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14752

GenBank Accession Number:

BC008304

GeneID (NCBI):

56902

UNIPROT ID:

Q9NRX1

Full Name:

partner of NOB1 homolog (S.
cerevisiae)

Calculated MW:

252 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

250134C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02143-1, PNO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85834-3-PBS. Detection antibody: 85834-2-PBS. Standard: Ag14752. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP02143-2, PNO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85834-2-PBS. Detection antibody: 85834-1-PBS. Standard: Ag14752. Range: 1.563-100 ng/mL.