

For Research Use Only

# PNO1 Recombinant antibody, PBS Only

Catalog Number: 85834-4-PBS



## Basic Information

Catalog Number:

85834-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14752

GenBank Accession Number:

BC008304

GeneID (NCBI):

56902

UNIPROT ID:

Q9NRX1

Full Name:

partner of NOB1 homolog (S. cerevisiae)

Calculated MW:

252 aa, 28 kDa

Observed MW:

28-35 kDa

Purification Method:

Protein A purification

CloneNo.:

250134F6

## Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

## Background Information

Human RNA binding gene partner of NOB1 (PNO1) was firstly isolated from human kidney and encodes 35-kDa protein, named PNO1. PNO1 as a ribosome biogenesis factor that contributes to 40S subunit assembly by binding with its partner NOB1. PNO1 is a highly conserved protein with a K homology (KH) domain at its C-terminal and two putative nuclear localization signals at its N-terminal. PNO1 also had great regulatory function in tumor metabolism. (PMID: 38071381, PMID: 33864661)

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85834-4-RR (PNO1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85834-4-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using PNO1 antibody (85834-4-RR, Clone: 250134F6 ) at dilution of 1:800 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85834-4-PBS in a different storage buffer formulation.