

For Research Use Only

# MBL Recombinant antibody

Catalog Number: 85835-4-RR



## Basic Information

Catalog Number:

85835-4-RR

Concentration:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3157

GenBank Accession Number:

NM\_010776.1

GeneID (NCBI):

17195

UNIPROT ID:

P41317

Full Name:

mannose-binding lectin (protein C) 2

Calculated MW:

26 kDa

Observed MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

250108E11

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : mouse liver tissue, rat liver tissue, mouse serum

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85835-4-RR (Mbl2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



rat liver tissue were subjected to SDS PAGE followed by western blot with 85835-4-RR (Mbl2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.