

For Research Use Only

Rat Fas Ligand/TNFSF6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85836-2-PBS



Basic Information

Catalog Number: 85836-2-PBS	GenBank Accession Number: NM_012908.1	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 25385	CloneNo.: 242892A9
Source: Rabbit	UNIPROT ID: P36940	
Isotype: IgG	Full Name: Fas ligand (TNF superfamily, member 6)	
Immunogen Catalog Number: EG3218	Calculated MW: 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2135-2, RAT Fas/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85836-2-PBS. Detection antibody: 85836-1-PBS. Standard: Eg3218. Range: 1.563-200 ng/mL.