

For Research Use Only

# Rat Fas Ligand/TNFSF6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85836-3-PBS



## Basic Information

<b>Catalog Number:</b> 85836-3-PBS	<b>GenBank Accession Number:</b> NM_012908.1	<b>Purification Method:</b> Protein A purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 25385	<b>CloneNo.:</b> 242892H8
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P36940	
<b>Isotype:</b> IgG	<b>Full Name:</b> Fas ligand (TNF superfamily, member 6)	
<b>Immunogen Catalog Number:</b> EG3218	<b>Calculated MW:</b> 31 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
rat

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MPO2135-1, RAT Fas/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85836-4-PBS. Detection antibody: 85836-3-PBS. Standard: Eg3218. Range: 1.563-200 ng/mL.