

For Research Use Only

SENP3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85837-2-PBS



Basic Information

Catalog Number:

85837-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11911

GenBank Accession Number:

BC080658

GeneID (NCBI):

26168

UNIPROT ID:

Q9H4L4

Full Name:

SUMO1/sentrin/SMT3 specific
peptidase 3

Calculated MW:

574 aa, 65 kDa

Purification Method:

Protein A purification

CloneNo.:

250153B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02141-2, SENP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85837-3-PBS. Detection antibody: 85837-2-PBS. Standard: Ag11911. Range: 0.781-100 ng/mL.