

For Research Use Only

SOX15 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85840-4-PBS



Basic Information

Catalog Number: 85840-4-PBS	GenBank Accession Number: BC000985	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 6665	CloneNo.: 250193D11
Source: Rabbit	UNIPROT ID: O60248	
Isotype: IgG	Full Name: SRY (sex determining region Y)-box 15	
Immunogen Catalog Number: AG22032	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02126-2, SOX15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85840-4-PBS. Detection antibody: 85840-3-PBS. Standard: Ag22032. Range: 0.469-60 ng/mL