

For Research Use Only

Mouse Tnfsf18 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85842-1-PBS



Basic Information

Catalog Number:

85842-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2850

GenBank Accession Number:

NM_183391.3

GeneID (NCBI):

240873

UNIPROT ID:

Q7TS55

Full Name:

tumor necrosis factor (ligand)
superfamily, member 18

Calculated MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

242937C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02134-2, MOUSE Tnfsf18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85842-3-PBS. Detection antibody: 85842-1-PBS. Standard: Eg2850. Range: 0.781-25 ng/mL.