

For Research Use Only

# VEGFR-1/FLT-1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 85844-3-PBS



## Basic Information

<b>Catalog Number:</b> 85844-3-PBS	<b>GenBank Accession Number:</b> NM_002019	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 2321	<b>CloneNo.:</b> 253459B11
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P17948-1	
<b>Immunogen Catalog Number:</b> EG7078	<b>Full Name:</b> fms-related tyrosine kinase 1 (vascular endothelial growth factor/vascular permeability factor receptor)	
	<b>Calculated MW:</b> 151kd	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Saliva of six individual healthy human donors was measured. The human VEGFR-1/FLT-1 concentration of detected samples was determined to be 4,798.9 ppg/mL with a range of 1,871.9-9,859.6 pg/mL.



Plasma of six individual healthy human donors was measured. The human VEGFR-1/FLT-1 concentration of detected samples was determined to be 571.1 pg/mL with a range of 405.2-721.5 pg/mL.



Sandwich ELISA standard curve of MP03426-1, Human VEGFR-1/FLT-1 Recombinant Matched Antibody Pair - PBS only. 85844-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg7078. 85844-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.