

For Research Use Only

ASF1B Recombinant antibody, PBS Only (Detector)

Catalog Number: 85848-2-PBS



Basic Information

Catalog Number:

85848-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17763

GenBank Accession Number:

BC010014

GeneID (NCBI):

55723

UNIPROT ID:

Q9NVP2

Full Name:

ASF1 anti-silencing function 1
homolog B (S. cerevisiae)

Calculated MW:

202 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

250195F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02155-1, ASF1B-specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85848-3-PBS. Detection antibody: 85848-2-PBS. Standard: Ag17763. Range: 0.781-100 ng/mL.