

For Research Use Only

DCUN1D1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85849-3-PBS



Basic Information

Catalog Number:

85849-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1440

GenBank Accession Number:

BC009478

GeneID (NCBI):

54165

UNIPROT ID:

Q96GG9

Full Name:

DCN1, defective in cullin neddylation
1, domain containing 1 (S. cerevisiae)

Calculated MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

250196E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02158-2, DCUN1D1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85849-3-PBS. Detection antibody: 85849-2-PBS. Standard: Ag1440. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02158-1, DCUN1D1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85849-3-PBS. Detection antibody: 85849-1-PBS. Standard: Ag1440. Range: 0.781-100 ng/mL.