

For Research Use Only

# PDRG1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85851-1-PBS



## Basic Information

|                                     |  |  |
|-------------------------------------|--|--|
| Catalog Number:<br>85851-1-PBS      | GenBank Accession Number:<br>BC009334        | Purification Method:<br>Protein A purification |
| Concentration:<br>1 mg/ml           | GeneID (NCBI):<br>81572                      | CloneNo.:<br>250203F8                          |
| Source:<br>Rabbit                   | UNIPROT ID:<br>Q9NUG6                        |  |
| Isotype:<br>IgG                     | Full Name:<br>p53 and DNA damage regulated 1 |  |
| Immunogen Catalog Number:<br>AG9109 | Calculated MW:<br>16 kDa                     |  |

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02151-2, PDRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85851-2-PBS. Detection antibody: 85851-1-PBS. Standard: Ag9109. Range: 0.781-100 ng/mL.