

For Research Use Only

# PDRG1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85851-2-PBS



## Basic Information

Catalog Number:

85851-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9109

GenBank Accession Number:

BC009334

GeneID (NCBI):

81572

UNIPROT ID:

Q9NUG6

Full Name:

p53 and DNA damage regulated 1

Calculated MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

250203G10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02151-1, PDRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85851-3-PBS. Detection antibody: 85851-2-PBS. Standard: Ag9109. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02151-2, PDRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85851-2-PBS. Detection antibody: 85851-1-PBS. Standard: Ag9109. Range: 0.781-100 ng/mL.