

For Research Use Only

PAEP/Glycodelin Recombinant antibody, PBS Only

Catalog Number: 85856-1-PBS



Basic Information

Catalog Number:

85856-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3037

GenBank Accession Number:

NM_002571.4

GeneID (NCBI):

5047

UNIPROT ID:

P09466-1

Full Name:

progesterone-associated endometrial protein

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

250206F4

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

