

For Research Use Only

# ACTN2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85857-2-PBS



## Basic Information

Catalog Number:

85857-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5459

GenBank Accession Number:

BC051770

GeneID (NCBI):

88

UNIPROT ID:

P35609

Full Name:

actinin, alpha 2

Calculated MW:

104 kDa

Purification Method:

Protein A purification

CloneNo.:

250229F7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02153-2, ACTN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85857-2-PBS. Detection antibody: 85857-1-PBS. Standard: Ag5459. Range: 0.781-50 ng/mL