

For Research Use Only

NARG1L Recombinant antibody, PBS Only (Detector)

Catalog Number: 85860-1-PBS



Basic Information

Catalog Number:

85860-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11492

GenBank Accession Number:

BC028112

GeneID (NCBI):

79612

UNIPROT ID:

Q6N069

Full Name:

NMDA receptor regulated 1-like

Calculated MW:

864 aa, 101 kDa

Purification Method:

Protein A purification

CloneNo.:

250181A1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02152-2, NARG1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85860-3-PBS. Detection antibody: 85860-1-PBS. Standard: Ag11492. Range: 0.781-100 ng/mL