

For Research Use Only

# NARG1L Recombinant antibody, PBS Only (Detector)

Catalog Number: 85860-2-PBS



## Basic Information

Catalog Number: 85860-2-PBS	GenBank Accession Number: BC028112	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 79612	CloneNo.: 250181E6
Source: Rabbit	UNIPROT ID: Q6N069	
Isotype: IgG	Full Name: NMDA receptor regulated 1-like	
Immunogen Catalog Number: AG11492	Calculated MW: 864 aa, 101 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02152-1, NARG1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85860-3-PBS. Detection antibody: 85860-2-PBS. Standard: Ag11492. Range: 0.781-100 ng/mL