

For Research Use Only

# NARG1L Recombinant antibody, PBS Only

Catalog Number: 85860-4-PBS



## Basic Information

Catalog Number:

85860-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11492

GenBank Accession Number:

BC028112

GeneID (NCBI):

79612

UNIPROT ID:

Q6N069

Full Name:

NMDA receptor regulated 1-like

Calculated MW:

864 aa, 101 kDa

Observed MW:

101 kDa

Purification Method:

Protein A purification

CloneNo.:

250181F4

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

NMDA receptor-regulated 1-like protein (NARG1L, also known as NAA16), is the auxiliary subunit of the N-terminal acetyltransferase A (NatA) complex which displays alpha (N-terminal) acetyltransferase activity.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85860-4-RR (NARG1L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85860-4-PBS in a different storage buffer formulation.