

For Research Use Only

Mouse Ccl19 Recombinant antibody, PBS Only (Detector)

Catalog Number:85861-2-PBS



Basic Information

Catalog Number: 85861-2-PBS	GenBank Accession Number: NM_011888.2	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 24047	CloneNo.: 242867C12
Source: Rabbit	UNIPROT ID: O70460	
Isotype: IgG	Full Name: chemokine (C-C motif) ligand 19	
Immunogen Catalog Number: EG3124	Calculated MW: 12 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02165-2, MOUSE CCL19/MIP-3 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85861-3-PBS. Detection antibody: 85861-2-PBS. Standard: Eg3124. Range: 0.781-100 ng/mL.