

For Research Use Only

Mouse B2m Recombinant antibody, PBS Only (Detector)

Catalog Number: 85871-1-PBS



Basic Information

Catalog Number:

85871-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009735.3

GeneID (NCBI):

12010

UNIPROT ID:

P01887

Full Name:

beta-2 microglobulin

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

250172D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02174-1, MOUSE Beta-2-Microglobulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85871-2-PBS. Detection antibody: 85871-1-PBS. Standard: Eg3086. Range: 0.313-20 ng/mL