

For Research Use Only

# Mouse B2m Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85871-2-PBS



## Basic Information

Catalog Number:

85871-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_009735.3

GeneID (NCBI):

12010

UNIPROT ID:

P01887

Full Name:

beta-2 microglobulin

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

250172G12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02174-2, MOUSE Beta-2-Microglobulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85871-3-PBS. Detection antibody: 85871-2-PBS. Standard: Eg3086. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP02174-1, MOUSE Beta-2-Microglobulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85871-2-PBS. Detection antibody: 85871-1-PBS. Standard: Eg3086. Range: 0.313-20 ng/mL