

For Research Use Only

Mouse Osm Recombinant antibody, PBS Only

Catalog Number:85873-3-PBS



Basic Information

Catalog Number: 85873-3-PBS	GenBank Accession Number: NM_001013365.2	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 18413	CloneNo.: 242891D11
Source: Rabbit	UNIPROT ID: P53347	
Isotype: IgG	Full Name: oncostatin M	
Immunogen Catalog Number: EG2921	Calculated MW: 30 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data