

For Research Use Only

CPLX1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85874-2-PBS



Basic Information

Catalog Number:

85874-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0387

GenBank Accession Number:

BC002471

GeneID (NCBI):

10815

UNIPROT ID:

O14810

Full Name:

complexin 1

Calculated MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

250131D1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02175-2, CPLX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85874-4-PBS. Detection antibody: 85874-2-PBS. Standard: Ag0387. Range: 6.25-200 ng/mL