

For Research Use Only

# GLUT5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85878-3-PBS



## Basic Information

Catalog Number:

85878-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26091

GenBank Accession Number:

BC001820

GeneID (NCBI):

6518

UNIPROT ID:

P22732

Full Name:

solute carrier family 2 (facilitated  
glucose/fructose transporter),  
member 5

Calculated MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

250150A6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02184-1, GLUT5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85878-3-PBS. Detection antibody: 85878-2-PBS. Standard: Ag26091. Range: 0.313-40 ng/mL