

For Research Use Only

Mouse Timp2 Recombinant antibody, PBS Only (Detector)

Catalog Number:85882-1-PBS



Basic Information

Catalog Number:

85882-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2960

GenBank Accession Number:

uniprot无NM号

GeneID (NCBI):

21858

UNIPROT ID:

P25785

Full Name:

tissue inhibitor of metalloproteinase

2

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

250231G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02176-1, MOUSE Timp2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85882-2-PBS. Detection antibody: 85882-1-PBS. Standard: Eg2960. Range: 0.781-100 ng/mL.