

For Research Use Only

# LGALS3BP Recombinant antibody, PBS Only (Detector)

Catalog Number: 85883-3-PBS



## Basic Information

Catalog Number:

85883-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3302

GenBank Accession Number:

NM\_005567

GeneID (NCBI):

3959

UNIPROT ID:

Q08380

Full Name:

lectin, galactoside-binding, soluble, 3 binding protein

Calculated MW:

585 aa, 65 kDa

Purification Method:

Protein A purification

CloneNo.:

250138B5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02172-2, LGALS3BP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85883-4-PBS. Detection antibody: 85883-3-PBS. Standard: Eg3302. Range: 0.781-50 ng/mL