

For Research Use Only

# RBBP6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85884-4-PBS



## Basic Information

Catalog Number:

85884-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2484

GenBank Accession Number:

BC029352

GeneID (NCBI):

5930

UNIPROT ID:

Q7Z6E9

Full Name:

retinoblastoma binding protein 6

Calculated MW:

1792 aa, 200 kDa

Purification Method:

Protein A purification

CloneNo.:

250142E4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02182-1, RBBP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85884-4-PBS. Detection antibody: 85884-3-PBS. Standard: Ag2484. Range: 0.625-80 ng/mL