

For Research Use Only

TAF1B Recombinant antibody, PBS Only (Capture)

Catalog Number: 85885-1-PBS



Basic Information

Catalog Number:

85885-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38211

GenBank Accession Number:

BC018137

GeneID (NCBI):

9014

UNIPROT ID:

Q53T94

Full Name:

TATA box binding protein (TBP)-
associated factor, RNA polymerase I,
B, 63kDa

Calculated MW:

588 aa, 69 kDa

Purification Method:

Protein A purification

CloneNo.:

250162C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02183-1, TAF1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85885-1-PBS. Detection antibody: 85885-2-PBS. Standard: Ag38211. Range: 0.781-100 ng/mL.