

For Research Use Only

WNK4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85886-3-PBS



Basic Information

Catalog Number:

85886-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17773

GenBank Accession Number:

BC136664

GeneID (NCBI):

65266

UNIPROT ID:

Q96J92

Full Name:

WNK lysine deficient protein kinase 4

Calculated MW:

1243 aa, 135 kDa

Purification Method:

Protein A purification

CloneNo.:

250156A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02181-1, WNK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85886-3-PBS. Detection antibody: 85886-2-PBS. Standard: Ag17773. Range: 0.469-30 ng/mL



Cytometric bead array standard curve of MP02181-2, WNK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85886-3-PBS. Detection antibody: 85886-1-PBS. Standard: Ag17773. Range: 0.469-60 ng/mL