

For Research Use Only

# Rat RANK/TNFRSF11A Recombinant antibody, PBS Only

Catalog Number: 85888-2-PBS



## Basic Information

Catalog Number:

85888-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3285

GenBank Accession Number:

NM\_001271235.1

GeneID (NCBI):

498206

UNIPROT ID:

F1M8Z6

Full Name:

Tumor necrosis factor receptor  
superfamily, member 11a

Calculated MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

242619D7

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

