

For Research Use Only

Rat RANK/TNFRSF11A Recombinant antibody, PBS Only (Capture)

Catalog Number: 85888-5-PBS



Basic Information

Catalog Number:

85888-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3285

GenBank Accession Number:

NM_001271235.1

GeneID (NCBI):

498206

UNIPROT ID:

F1M8Z6

Full Name:

Tumor necrosis factor receptor
superfamily, member 11a

Calculated MW:

68 kDa

Purification Method:

Protein A purification

CloneNo.:

242619F9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02178-3, Human TNFRSF11A Recombinant Matched Antibody Pair - PBS only. 85888-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3285. 85888-1-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean TNFRSF11A concentration was determined to be 15.90 ng/mL in rat spleen extract based on a 2.70 mg/mL extract load.