

For Research Use Only

ACADVL Recombinant antibody

Catalog Number: 85889-1-RR



Basic Information

Catalog Number:

85889-1-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6016

GenBank Accession Number:

BC012912

GeneID (NCBI):

37

UNIPROT ID:

P49748

Full Name:

acyl-Coenzyme A dehydrogenase,
very long chain

Calculated MW:

70 kDa

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

250199C12

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, HepG2 cells, HEK-293T cells, A431
cells, L02 cells, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85889-1-RR (ACADVL antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.