

For Research Use Only

ASNS Recombinant antibody, PBS Only (Detector)

Catalog Number: 85892-3-PBS



Basic Information

Catalog Number: 85892-3-PBS	GenBank Accession Number: BC014621	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 440	CloneNo.: 250177E8
Source: Rabbit	UNIPROT ID: P08243	
Isotype: IgG	Full Name: asparagine synthetase	
Immunogen Catalog Number: AG6393	Calculated MW: 64 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02189-2, ASNS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85892-4-PBS. Detection antibody: 85892-3-PBS. Standard: Ag6393. Range: 0.938-120 ng/mL.