

For Research Use Only

ASNS Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85892-4-PBS



Basic Information

Catalog Number:

85892-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6393

GenBank Accession Number:

BC014621

GeneID (NCBI):

440

UNIPROT ID:

P08243

Full Name:

asparagine synthetase

Calculated MW:

64 kDa

Purification Method:

Protein A purification

CloneNo.:

250177H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02189-2, ASNS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85892-4-PBS. Detection antibody: 85892-3-PBS. Standard: Ag6393. Range: 0.938-120 ng/mL.



Cytometric bead array standard curve of MP02189-1, ASNS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85892-2-PBS. Detection antibody: 85892-4-PBS. Standard: Ag6393. Range: 0.938-120 ng/mL.