

For Research Use Only

# TGFBR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85899-3-PBS



## Basic Information

Catalog Number:

85899-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2943

GenBank Accession Number:

NM\_003242.5

GeneID (NCBI):

7048

UNIPROT ID:

P37173-1

Full Name:

transforming growth factor, beta receptor II (70/80kDa)

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

250238E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02190-2, TGFBR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85899-3-PBS. Detection antibody: 85899-4-PBS. Standard: Eg2943. Range: 0.313-40 ng/mL