

For Research Use Only

TGFBR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85899-4-PBS



Basic Information

Catalog Number:

85899-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2943

GenBank Accession Number:

NM_003242.5

GeneID (NCBI):

7048

UNIPROT ID:

P37173-1

Full Name:

transforming growth factor, beta receptor II (70/80kDa)

Calculated MW:

65 kDa

Purification Method:

Protein A purification

CloneNo.:

250238A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02190-2, TGFBR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85899-3-PBS. Detection antibody: 85899-4-PBS. Standard: Eg2943. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP02190-1, TGFBR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85899-2-PBS. Detection antibody: 98438-1-PBS. Standard: Eg2943. Range: 0.313-20 ng/mL