

For Research Use Only

SATB2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85922-7-PBS



Basic Information

Catalog Number:

85922-7-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_015265

GeneID (NCBI):

23314

UNIPROT ID:

Q9UPW6

Full Name:

SATB homeobox 2

Calculated MW:

83 kDa

Purification Method:

Protein A purification

CloneNo.:

260789A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean SATB2 concentration was determined to be 2.3 ng/mL in THP-1 cell extract based on a 1.0 mg/mL extract load, 0.6 ng/mL in K-562 cell extract based on a 1.0 mg/mL extract load, 0.1 ng/mL in A2780 cell extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03557-1, SATB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85922-7-PBS. Detection antibody: 85922-6-PBS. Standard: Ag34010. Range: 0.156-10 ng/mL.