

For Research Use Only

# TH Recombinant antibody, PBS Only (Detector)

Catalog Number: 85929-2-PBS



## Basic Information

Catalog Number: 85929-2-PBS	GenBank Accession Number: BC104967	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 7054	CloneNo.: 250146C6
Source: Rabbit	UNIPROT ID: P07101	
Isotype: IgG	Full Name: tyrosine hydroxylase	
Immunogen Catalog Number: AG23034	Calculated MW: 528 aa, 59 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean TH concentration was determined to be 107.4 ng/mL in Mouse brain tissue extract based on a 7.0 mg/mL extract load.



Cytometric bead array standard curve of MP02201-1, TH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85929-3-PBS. Detection antibody: 85929-2-PBS. Standard: Ag23034. Range: 3.125-200 ng/mL.