

For Research Use Only

TH Recombinant antibody, PBS Only (Capture)

Catalog Number: 85929-3-PBS



Basic Information

Catalog Number: 85929-3-PBS	GenBank Accession Number: BC104967	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 7054	CloneNo.: 250146A3
Source: Rabbit	UNIPROT ID: P07101	
Isotype: IgG	Full Name: tyrosine hydroxylase	
Immunogen Catalog Number: AG23034	Calculated MW: 528 aa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TH concentration was determined to be 107.4 ng/mL in Mouse brain tissue extract based on a 7.0 mg/mL extract load.



Cytometric bead array standard curve of MP02201-1, TH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85929-3-PBS. Detection antibody: 85929-2-PBS. Standard: Ag23034. Range: 3.125-200 ng/mL.