

For Research Use Only

# FLI1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85931-2-PBS



## Basic Information

Catalog Number:

85931-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36169

GenBank Accession Number:

BC010115

GeneID (NCBI):

2313

UNIPROT ID:

Q01543

Full Name:

Friend leukemia virus integration 1

Calculated MW:

452 aa, 51 kDa

Purification Method:

Protein A purification

CloneNo.:

250294A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean FLI1 concentration was determined to be 334.3 ng/mL in Jurkat cell extract based on a 1.9 mg/mL extract load.



Cytometric bead array standard curve of MP02198-2, FLI1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85931-3-PBS. Detection antibody: 85931-2-PBS. Standard: Ag36169. Range: 1.563-100 ng/mL.