

For Research Use Only

# DGUOK Recombinant antibody

Catalog Number: 85932-2-RR



## Basic Information

Catalog Number:

85932-2-RR

Concentration:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9361

GenBank Accession Number:

BC015757

GeneID (NCBI):

1716

UNIPROT ID:

Q16854

Full Name:

deoxyguanosine kinase

Calculated MW:

277 aa, 32 kDa

Observed MW:

28-32 kDa

Purification Method:

Protein A purification

CloneNo.:

250161D11

Recommended Dilutions:

WB: 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, A549 cells, U-87 MG cells, K-562 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85932-2-RR (DGUOK antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.