

For Research Use Only

NUMB Recombinant antibody, PBS Only (Capture)

Catalog Number: 85939-3-PBS



Basic Information

Catalog Number:

85939-3-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31553

GenBank Accession Number:

NM_001005743

GeneID (NCBI):

8650

UNIPROT ID:

P49757

Full Name:

numb homolog (Drosophila)

Calculated MW:

71KD

Purification Method:

Protein A purification

CloneNo.:

250417D8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean NUMB concentration was determined to be 8.7 ng/mL in A431 cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02209-1, NUMB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85939-3-PBS. Detection antibody: 85939-2-PBS. Standard: Ag31553. Range: 0.625-80 ng/mL.



Biolayer Interferometry (BLI) kinetic assay of 85939-3-PBS against Human NUMB was performed. The affinity constant is 0.738 nM.