

For Research Use Only

NUMB Recombinant antibody, PBS Only (Detector)

Catalog Number: 85939-5-PBS



Basic Information

Catalog Number:

85939-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31553

GenBank Accession Number:

NM_001005743

GeneID (NCBI):

8650

UNIPROT ID:

P49757

Full Name:

numb homolog (Drosophila)

Calculated MW:

71KD

Purification Method:

Protein A purification

CloneNo.:

250417D2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02209-2, Human NUMB Recombinant Matched Antibody Pair - PBS only. 85939-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag31553. 85939-5-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean NUMB concentration was determined to be 20.70 ng/mL in HeLa cell extract based on a 1.30 mg/mL extract load.