

For Research Use Only

Mouse Icam5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85941-2-PBS



Basic Information

Catalog Number:

85941-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2868

GenBank Accession Number:

N/A

GeneID (NCBI):

15898

UNIPROT ID:

Q60625.1

Full Name:

intercellular adhesion molecule 5,
telencephalin

Calculated MW:

97 kDa

Purification Method:

Protein A purification

CloneNo.:

243100D8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ICAM5 concentration was determined to be 630.0 ng/mL in Mouse brain tissue extract based on a 7 mg/mL extract load.



Cytometric bead array standard curve of MP02210-1, MOUSE ICAM5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85941-3-PBS. Detection antibody: 85941-2-PBS. Standard: Eg2868. Range: 1.563-200 ng/mL.