

For Research Use Only

IL-23 Recombinant antibody

Catalog Number: 85947-2-RR



Basic Information

Catalog Number:

85947-2-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1254

GenBank Accession Number:

GeneID (NCBI):

51561

UNIPROT ID:

Q9NPF7

Full Name:

interleukin 23, alpha subunit p19

Purification Method:

Protein A purification

CloneNo.:

250414G4

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Recombinant protein,

Background Information

Storage

Storage:

Store at -20°C . Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85947-2-RR (IL-23 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Recombinant protein were subjected to SDS PAGE followed by western blot with 85947-2-RR (IL-23 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.