

For Research Use Only

MCEE Recombinant antibody, PBS Only (Detector)

Catalog Number: 85948-2-PBS



Basic Information

Catalog Number:

85948-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9665

GenBank Accession Number:

BC020825

GeneID (NCBI):

84693

UNIPROT ID:

Q96PE7

Full Name:

methylmalonyl CoA epimerase

Calculated MW:

176 aa, 19 kDa

Purification Method:

Protein A purification

CloneNo.:

250393E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MCEE concentration was determined to be 644.5 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02207-1, MCEE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85948-3-PBS. Detection antibody: 85948-2-PBS. Standard: Ag9665. Range: 0.781-100 ng/mL.