

For Research Use Only

HOXA10 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:85951-3-PBS



Basic Information

Catalog Number: 85951-3-PBS	GenBank Accession Number: BC013971	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3206	CloneNo.: 250418C7
Isotype: IgG	UNIPROT ID: P31260	
Immunogen Catalog Number: AG24156	Full Name: homeobox A10	
	Calculated MW: 393 aa, 41 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HOXA10 concentration was determined to be 21.8 ng/mL in K-562 cell extract based on a 3.0 mg/mL extract load and 56.2 ng/mL in HEK-293T cell extract based on a 1.3 mg/mL extract load.



Sandwich ELISA standard curve of MP02388-1, Human HOXA10 Recombinant Matched Antibody Pair - PBS only. 85951-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag24156. 85951-3-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL