

For Research Use Only

# ANXA13 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 85955-2-PBS



## Basic Information

Catalog Number:

85955-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18581

GenBank Accession Number:

BC125158

GeneID (NCBI):

312

UNIPROT ID:

P27216

Full Name:

annexin A13

Calculated MW:

316 aa, 35 kDa

Purification Method:

Protein A purification

CloneNo.:

250350D3

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Annexin A13 concentration was determined to be 46.6 ng/mL in Jurkat cell extract based on a 1.3 mg/mL extract load.



Sandwich ELISA standard curve of MP02684-1, Human Annexin A13 Recombinant Matched Antibody Pair - PBS only. 85955-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag18581. 85955-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.